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Table 1. shows the summary of the physical properties of some typical cryogens. These possess low boiling point, relatively lower latent heat, large volume ratio of the gas and liquid. In using the cryogen, there exist possibilities of accidents of frostbite, lack of oxygen (suffocation), explosion, and others. It is important to understand their properties very well and handle the cryogens.

Table 1. 寒剤の性質 Cryogen properties

液体 Liquid	分子量 Molecular weight	沸点 (K) Boiling point	密度 (g/cm ³) [†] Density	潜熱 (cal/cm ³) Latent heat	容積比* Volume ratio
O ₂	32.0	90.19	1.144	58.2	800
N ₂	28.01	77.35	0.808	38.5	647
⁴ He	4.003	4.216	0.125	0.62	700
³ He	3.000	3.1905	0.059	0.12	440

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Volume ratio: the ratio between the volume of the liquid at the boiling point and that of the gas at 0°C and 1 atm.

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In the handling the cryogen, your hands and foots may be frostbitten by causes mentioned below.

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When the bare hands or body are exposed to a cryogen or cold gas.

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When you touch a cold metal such as the surface on a delivery tube. When a small amount of liq. N₂ drops on your hands, your hands protected by the N₂ gas is not frostbitten. However, when a cold metal cooled by the liq. N₂ is touched by bare hands, the hands and metal sticks together

because of the no protection of the gas, and then, the hands are heavily frostbitten. In the cold metal cooled by the liq. He, there becomes a serious situation.(See reference materials p.22)

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When you touch cryogenics directly by your wet hands or feet.

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You should use **leather gloves** or **thermal insulating gloves** in the handling the cryogen. Do not use cotton or woolen work gloves in which liquid N₂ sinks.

Fig. 2. 断熱手袋，革製手袋と軍手 Thermal insulating gloves, leather gloves and cotton work gloves

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Please brush away the Liq. N₂ quickly when the liq. N₂ splashes.

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You should use goggles to prevent the liq. N₂ from getting into the eye when the Liq. N₂ splashes.

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You should soak the frostbitten skin in a hot water. Do not use a drier. After that, go to a hospital.

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When your hand stick to a cold metal, you warm up the metal up to the temperature of tepid water, and them, take your hand off the metal slowly. After that, you should go to a hospital.

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When the Liq. N₂ enters your eyes, you should wash your eyes by clean water.

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Vaporization of liquid N₂ will increase its volume and decrease the concentration of O₂ in the room.

When the cryogenics are used in the closed room, we must ventilate the room.

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The rapid oxygen-shortage is caused by breathing a gas without oxygen. **When you come in the room in which the oxygen is 0%, you fall down in one breathing, and have no chance of resuscitation.** To breathe near the evaporating gas of liq. N₂ is a suicidal act.

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The gradual oxygen-shortage is caused slowly by breathing an air in which the oxygen concentration decreases gradually. (See Table 2.)

Table 2. 酸素濃度低下(緩急性酸欠)による人体への影響
Effects on the Human Body due to Gradual Oxygen-shortage

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[酸欠事故の被害拡大の要因] Factors of the Accident in Oxygen-shortage

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Secondary Disaster: There is a high possibility that a rescuer coming into the room without oxygen suffer damage because he can not find why a person falls down in the room without smoke or irritating smell.

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Lack of Judgment Ability: A person coming into the room without oxygen can not understand that his headache, nausea, and coma are attributed to the lack of oxygen.

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Sharp Drop in Resuscitation Opportunity: A chance of resuscitation decreases rapidly in 3~5 minutes after the stop of breathing, and becomes zero above 7 minutes (Fig. 3). Even if a person is resuscitated, he has sometimes after effects of disease such as speech impediment, motility disturbance, narrowing of visual field, hallucination, neurosis.

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A rapid treatment preventing the secondary disaster is required.

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When you use liq. N₂ in a closed room, you should release the N₂ gas out the room by using a ventilating fan or exhaust pipe. (See reference materials p.21)

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When you take a vessel of liq. N₂ by car, you must open the windows.

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When liq. N₂ or He is transported by elevator, the vessels should travel alone in the elevator, without any people accompanying them inside the elevator. (See reference materials p.19, 20)

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When the inside of a room seems to be in the lack of oxygen, you should ventilate the room with a fresh air, and you may enter the room after you find that the oxygen concentration becomes 20% or more.

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Liq. He and liq. N₂ boil even in the vessels because the heat enters from outside. When the liq. He and N₂ are transformed into gases at 273 K and 1 atm, their volumes expand up to 700 and 647 times, respectively. When cryogen vessels are sealed, there is a danger of explosion because the pressure in the vessels becomes very high. **Do not seal cryogen vessels. (See Table 1.)**

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The vent of the vessel containing a cryogen is closed by ice.

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All of valves in the sealed vessel are closed.

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You should mount check valve in the vent by preventing air entering into the vessels.

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When you use relatively large sealed vessels with some relief valves, please read those manuals and understand how to use them very well.

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—The Neck in the Vessel is Weak Mechanically

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The neck in the vessel is weak mechanically. Therefore, do not knock the vessel down, and put the large vibration.

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—Notes for Using the Vessels and Transfer of the Cryogenics


Fig. 5. 液体窒素シールド型ヘリウム容器
Liq. He vessel shielded by liq. N₂

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If there are some contaminations and water in the vessel, we should remove them before the transfer of cryogens.

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It is not good that cryogen is poured suddenly into a vessel at room temperature.

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There is a possibility that the large vessel over 100 L stumbles over a stone.

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Please open the cap of the closed vessel after the inner pressure is released.

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Do not use a safety valve as a release valve.

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Do not leave the place where we transfer cryogens. In the transfer of cryogens, we have possibility to be in unexpected accidents. We always need to watch circumstances carefully.

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Do not handle anything except transfer valve in liquid N₂ tank. There are some relief and safety valves in liquid N₂ tank. If a user handles these valves, there is possibility to increase the pressure in the tank. **Please handle only transfer valve in liquid N₂ tank.** If you find some anomalies, please contact the staff of Analysis Center of Life Science.

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Close the valve tightly after the transfer of cryogenics. The valve is frozen because the cryogenics come out little by little.

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—Pay Attention to Safety and Suppression of Cryogen Losses in the Transfer of Cryogenics

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In the transfer of cryogenics, we should carry out precooling by the cold gas. At first, we put the transfer tube into the bottom of the vessel, and cool slowly the tube and inner part of the vessel.

寒 剤 Cryogenics	He		N ₂
冷却温度 Cool down temperature (K)	4.2		77.3
初めの温度 Initial temperature (K)	300	77	300
潜熱だけを利用 Case for using the latent heat of the gas	31	2.2	0.46
潜熱とガスエンタルピーを利用 Case for using the latent heat and enthalpy of the gas	0.79	0.15	0.29

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When liquid N_2 is exposed to the air, the oxygen is condensed in the liquid N_2 . We have the hazard of explosion when flammable gases or organic compounds come into the liquid N_2 with the oxygen. In order to avoid a back current of the air, please use Bunsen valve or bubbler.

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Notes for transporting the containers of Liq. N₂ by elevator

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When liquid nitrogen is transported by elevator, the containers should travel alone in the elevator, without any people accompanying them inside the elevator. Also, the containers should bear a warning notice of the prescribed format to prevent anyone from getting in the elevator with them.

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When liquid nitrogen is transported by hand in a small portable container, the person carrying should use stairs wherever possible. If an elevator has to be used, the container should bear a notice of the prescribed format on its front so that nobody will get in the elevator with it.

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
Note 1) The prescribed format for the notice can be downloaded from the web page of the Natural Science Center for Basic Research and Development's Low Temperature Instrumental Analysis Department at <http://www.sci.hiroshima-u.ac.jp/kiki/>, or from the web page of the Safety and Health Committee at <http://home.hiroshima-u.ac.jp/anzen/>.

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